

兔抗 SLC27A2 多克隆抗体

中文名称：兔抗 SLC27A2 多克隆抗体

英文名称： Anti-SLC27A2 rabbit polyclonal antibody

别名： VLCS; FATP2; VLACS; ACSVL1; FACVL1; hFACVL1; HsT17226

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： SLC27A2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格 Background:	The protein encoded by this gene is an isozyme of long-chain fatty-acid-coenzyme A ligase family. Although differing in substrate specificity, subcellular localization, and tissue distribution, all isozymes of this family convert free long-chain fatty acids into fatty acyl-CoA esters, and thereby play a key role in lipid biosynthesis and fatty acid degradation. This isozyme activates long-chain, branched-chain and very-long-chain fatty acids containing 22 or more carbons to their CoA derivatives.
Applications:	ELISA, IHC
Name of antibody:	SLC27A2
Immunogen:	Fusion protein of human SLC27A2
Full name:	solute carrier family 27 (fatty acid transporter), member 2

Synonyms :	VLCS; FATP2; VLACS; ACSVL1; FACVL1; hFACVL1; Hs T17226
SwissProt:	O14975
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

