

兔抗 ZG16 多克隆抗体

中文名称：兔抗 ZG16 多克隆抗体

英文名称： Anti-ZG16 rabbit polyclonal antibody

别名： zymogen granule protein 16; JCLN; JCLN1; ZG16A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ZG16

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	May play a role in protein trafficking. May act as a linker molecule between the submembranous matrix on the luminal side of zymogen granule membrane (ZGM) and aggregated secretory proteins during granule formation in the TGN.
Applications:	ELISA, WB
Name of antibody:	ZG16
Immunogen:	Fusion protein of human ZG16
Full name:	zymogen granule protein 16
Synonyms:	JCLN; JCLN1; ZG16A
SwissProt:	O60844

ELISA Recommended dilution:	5000-10000
WB Predicted band size:	18 kDa
WB Positive control:	Human sigmoid tissue, Mouse large intestine tissue, Mouse small intestines tissue, Mouse Pancreas tissue lysates
WB Recommended dilution:	500-2000

