

兔抗 KLK11 多克隆抗体

中文名称:	兔抗 KLK11 多克隆抗体
英文名称:	Anti-KLK11 rabbit polyclonal antibody
别名:	TLSP; PRSS20
相关类别:	一抗
抗原:	KLK11
储存:	冷冻 (-20℃)
宿主:	Rabbit
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Alternate splicing of this gene results in multiple transcript variants encoding distinct isoforms which are differentially expressed.
Applications:	ELISA, IHC
Name of antibody:	KLK11
Immunogen:	Fusion protein of human KLK11

Full name:	kallikrein-related peptidase 11
Synonyms :	TLSP; PRSS20
SwissProt:	Q9UBX7
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

