

## 兔抗 POTE<sub>G</sub>/H 多克隆抗体

中文名称: 兔抗 POTE<sub>G</sub>/H 多克隆抗体

英文名称: Anti-POTE<sub>G</sub>/H rabbit polyclonal antibody

别名: A26C2; ACTBL1; POTE14; POTE22; CT104.4; POTE-14; POTE14alpha/A26C3; ACTBL1; POTE22; CT104.7

抗原: POTE<sub>G</sub>/H

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Belongs to the POTE family. 3 isoforms of the human protein are produced by alternative splicing.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	POTE <sub>G</sub> /H
<b>Immunogen:</b>	Synthesized peptide derived from internal of human A26C2/3.
<b>Full name:</b>	POTE ankyrin domain family, member G/H
<b>Synonyms :</b>	A26C2; ACTBL1; POTE14; POTE22; CT104.4; POTE-14; POTE14alpha/A26C3; ACTBL1; POTE22; CT104.7
<b>SwissProt:</b>	Q6S5H5/Q6S545

<b>IHC positive control:</b>	Human prostate carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000



