

兔抗 KLK12 多克隆抗体

中文名称：兔抗 KLK12 多克隆抗体

英文名称：Anti-KLK12 rabbit polyclonal antibody

级别：KLKL5; KLK-L5

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：KLK12

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

ackground:	Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Alternate splicing of this gene results in three transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	KLK12
Immunogen:	Synthetic peptide of human KLK12

Full name:	kallikrein related peptidase 12
Synonyms:	KLKL5; KLK-L5
SwissProt:	Q9UKR0
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	10-50

kDa

95 —

72 —

55 —

36 —

28 —

17 —

10 —



