

## 兔抗 ZGPAT 多克隆抗体

中文名称：兔抗 ZGPAT 多克隆抗体

英文名称：Anti-ZGPAT rabbit polyclonal antibody

别名：zinc finger CCCH-type and G-patch domain containing; ZIP; ZC3H9; GPATC6; GPATCH6; ZC3HDC9; KIAA1847

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZGPAT

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcription repressor that specifically binds the 5' -GGAG[GA]A[GA]A-3' consensus sequence. Represses transcription by recruiting the chromatin multiprotein complex NuRD to target promoters. Negatively regulates expression of EGFR, a gene involved in cell proliferation, survival and migration. Its ability to repress genes of the EGFR pathway suggest it may act as a tumor suppressor. Able to suppress breast carcinogenesis.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	ZGPAT
<b>Immunogen:</b>	Fusion protein of human ZGPAT
<b>Full name:</b>	zinc finger CCCH-type and G-patch domain containing
<b>Synonyms:</b>	ZIP; ZC3H9; GPATC6; GPATCH6; ZC3HDC9; KIAA1847
<b>SwissProt:</b>	Q8N5A5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	40-200





