

## 兔抗 SLC29A4 多克隆抗体

中文名称：兔抗 SLC29A4 多克隆抗体

英文名称： Anti-SLC29A4 rabbit polyclonal antibody

别名： ENT4; PMAT

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC29A4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a member of the SLC29A/ENT transporter protein family. The encoded membrane protein catalyzes the reuptake of monoamines into presynaptic neurons, thus determining the intensity and duration of monoamine neural signaling. It has been shown to transport several compounds, including serotonin, dopamine, and the neurotoxin 1-methyl-4-phenylpyridinium. Alternative splicing results in multiple transcript variants.

<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SLC29A4
<b>Immunogen:</b>	Fusion protein of human SLC29A4
<b>Full name:</b>	solute carrier family 29 (equilibrative nucleoside transporter), member 4
<b>Synonyms :</b>	ENT4; PMAT
<b>SwissProt:</b>	Q7RTT9
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	293T cells
<b>WB Recommended dilution:</b>	500-2000

