

## 兔抗 POU4F1 多克隆抗体

中文名称：兔抗 POU4F1 多克隆抗体

英文名称：Anti-POU4F1 rabbit polyclonal antibody

别名：BRN3A; RDC-1; Oct-T1; brn-3A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：POU4F1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the POU-IV class of neural transcription factors. This protein is expressed in a subset of retinal ganglion cells and may be involved in the developing sensory nervous system. This protein may also promote the growth of cervical tumors. A translocation of this gene is associated with some adult acute myeloid leukemias.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	POU4F1
<b>Immunogen:</b>	Synthetic peptide of human POU4F1
<b>Full name:</b>	POU class 4 homeobox 1p

<b>Synonyms :</b>	BRN3A; RDC-1; Oct-T1; brn-3A
<b>SwissProt:</b>	Q01851
<b>ELISA Recommended dilution :</b>	1000-5000
<b>IHC positive control:</b>	Human Lymphoma
<b>IHC Recommend dilution:</b>	25-100

