

兔抗 POU4F3 多克隆抗体

中文名称：兔抗 POU4F3 多克隆抗体

英文名称： Anti-POU4F3 rabbit polyclonal antibody

别名： BRN3C; DFNA15

抗原： POU4F3

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the POU-domain family of transcription factors. POU-domain proteins have been observed to play important roles in control of cell identity in several systems. This protein is found in the retina and may play a role in determining or maintaining the identities of a small subset of visual system neurons. Defects in this gene are the cause of non-syndromic sensorineural deafness autosomal dominant type 15.
Applications:	WB
Name of antibody:	POU4F3
Immunogen:	Synthesized peptide derived from C-terminal of human PO

	U4F3.
Full name:	POU class 4 homeobox 3
Synonyms :	BRN3C; DFNA15
SwissProt:	Q15319
WB Predicted band size:	37 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

