

## 兔抗 SLC32A1 多克隆抗体

中文名称：兔抗 SLC32A1 多克隆抗体

英文名称：Anti-SLC32A1 rabbit polyclonal antibody

别名：VGAT; VIAAT

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC32A1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is an integral membrane protein involved in gamma-aminobutyric acid (GABA) and glycine uptake into synaptic vesicles. The encoded protein is a member of amino acid/polyamine transporter family II.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC32A1
<b>Immunogen:</b>	Synthetic peptide of human SLC32A1
<b>Full name:</b>	solute carrier family 32 (GABA vesicular transporter), member 1
<b>Synonyms :</b>	VGAT; VIAAT

<b>SwissProt:</b>	Q9H598
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

