

兔抗 ZKSCAN5 多克隆抗体

中文名称：兔抗 ZKSCAN5 多克隆抗体

英文名称： Anti-ZKSCAN5 rabbit polyclonal antibody

别名： ZFP95; ZNF914; ZSCAN37

抗原： ZKSCAN5

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a zinc finger protein of the Kruppel family. The protein contains a SCAN box and a KRAB A domain. A similar protein in mouse is differentially expressed in spermatogenesis. Two alternatively spliced transcript variants differing only in the 5' UTR have been described. Additional variants have been found, but their full-length sequences have not been determined.
Applications:	IF
Name of antibody:	ZKSCAN5
Immunogen:	Synthesized peptide derived from internal of human ZNF95.

Full name:	zinc finger with KRAB and SCAN domains 5
Synonyms :	ZFP95; ZNF914; ZSCAN37
SwissProt:	Q9Y2L8
IF positive control:	A549 cells
IF Recommend dilution:	100-500

