

兔抗 PPA1 多克隆抗体

中文名称：兔抗 PPA1 多克隆抗体

英文名称：Anti-PPA1 rabbit polyclonal antibody

别名：IOPPP; MGC111556; PP; PP1; SID6-8061

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：PPA1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	PPA1(Inorganic pyrophosphatase) is also named as IOPPP, PP and PPase. It maintains the thermodynamic favorability of these reactions by catalyzing the hydrolysis of pyrophosphates into organic phosphates, which are then exported across the cell membrane. The erythrocyte inorganic pyrophosphatase molecule would appear to be a dimer.
Applications:	WB
Name of antibody:	PPA1
Immunogen:	Fusion protein of human PPA1
Full name:	pyrophosphatase (inorganic) 1
Synonyms :	IOPPP; MGC111556; PP; PP1; SID6-8061

SwissProt:	Q15181
WB Predicted band size:	33 kDa
WB Positive control:	Human heart tissue
WB Recommended dilution:	500-2000

