

兔抗 SLC35E2A 多克隆抗体

中文名称: 兔抗 SLC35E2A 多克隆抗体

英文名称: Anti-SLC35E2A rabbit polyclonal antibody

别名: solute carrier family 35 member E2A; SLC35E2

相关类别: 一抗

抗原: SLC35E2A

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Putative transporter.
Applications:	ELISA, IHC
Name of antibody:	SLC35E2A
Immunogen:	Fusion protein of human SLC35E2A
Full name:	solute carrier family 35 member E2A
Synonyms:	SLC35E2
SwissProt:	P0CK97
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	25-100

