

兔抗 ZMAT2 多克隆抗体

中文名称：兔抗 ZMAT2 多克隆抗体

英文名称：Anti-ZMAT2 rabbit polyclonal antibody

别名：Snu23; Ptg-12; hSNU23

相关类别：一抗

抗原：ZMAT2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugat

技术规格

Background:	The function of ZMAT2 remains unknown.
Applications:	ELISA, WB, IHC
Name of antibody:	ZMAT2
Immunogen:	Full length fusion protein
Full name:	zinc finger, matrin-type 2
Synonyms :	Snu23; Ptg-12; hSNU23
SwissProt:	Q96NC0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and human breast cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	24 kDa

WB Positive control:	Jurkat,231,Hela,k562 and A172 cell,Mouse brain tissue lysates
WB Recommended dilution:	200-1000



