

## 兔抗 SLC35F6 多克隆抗体

中文名称：兔抗 SLC35F6 多克隆抗体

英文名称： Anti-SLC35F6 rabbit polyclonal antibody

别名： solute carrier family 35 member F6; ANT2BP; TANGO9; C2orf18

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC35F6

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | Involved in the maintenance of mitochondrial membrane potential in pancreatic ductal adenocarcinoma (PDAC) cell s. Promotes pancreatic ductal adenocarcinoma (PDAC) cell growth. May play a role as a nucleotide-sugar transporter . |
| <b>Applications:</b>     | ELISA, IHC   |
| <b>Name of antibody:</b> | SLC35F6  |
| <b>Immunogen:</b>        | Fusion protein of human SLC35F6  |
| <b>Full name:</b>        | solute carrier family 35 member F6   |
| <b>Synonyms:</b>         | ANT2BP; TANGO9; C2orf18  |
| <b>SwissProt:</b>        | Q8N357   |

|                                    |  |
|------------------------------------|--|
| <b>ELISA Recommended dilution:</b> | 5000-10000                               |
| <b>IHC positive control:</b>       | Human colorectal cancer and Human tonsil |
| <b>IHC Recommend dilution:</b>     | 100-300                                  |



