

兔抗 PPAT 多克隆抗体

中文名称：兔抗 PPAT 多克隆抗体

英文名称： Anti-PPAT rabbit polyclonal antibody

别名： GPAT; PRAT; ATASE

相关类别： 一抗

抗原： PPAT

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the purine/pyrimidine phosphoribosyltransferase family. It is a regulatory allosteric enzyme that catalyzes the first step of de novo purine nucleotide biosynthetic pathway. This gene and PAICS/AIRC gene, a bifunctional enzyme catalyzing steps six and seven of this pathway, are located in close proximity on chromosome 4, and divergently transcribed from an intergenic region.
Applications:	ELISA, WB, IHC
Name of antibody:	PPAT
Immunogen:	Fusion protein of human PPAT

Full name:	phosphoribosyl pyrophosphate amidotransferase
Synonyms :	GPAT; PRAT; ATASE
SwissProt:	Q06203
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and Human brain
IHC Recommend dilution:	100-300
WB Predicted band size:	57 kDa
WB Positive control:	Hela, Raji, Jurkat, A549 and NIH/3T3 cells
WB Recommended dilution:	500-2000



