

兔抗 KNOP1 多克隆抗体

中文名称: 兔抗 KNOP1 多克隆抗体

英文名称: Anti-KNOP1 rabbit polyclonal antibody

别名: TSG118; FAM191A; 101F10.1; C16orf88

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

克隆类型: rabbit polyclonal

抗原: KNOP1

反应种属: Human, Mouse

标记物: Unconjugate

技术规格

Background:	poly(A) RNA binding
Applications:	ELISA, IHC
Name of antibody:	KNOP1
Immunogen:	Synthetic peptide of human KNOP1
Full name:	lysine-rich nucleolar protein 1
Synonyms:	TSG118; FAM191A; 101F10.1; C16orf88
SwissProt:	Q1ED39

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human liver cancer
IHC Recommend dilution:	40-200



