

兔抗 PPCDC 多克隆抗体

中文名称：兔抗 PPCDC 多克隆抗体

英文名称：Anti-PPCDC rabbit polyclonal antibody

别名：phosphopantothenoylcysteine decarboxylase; coaC; MDS018; PPC-DC

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PPCDC

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Biosynthesis of coenzyme A (CoA) from pantothenic acid (vitamin B5) is an essential universal pathway in prokaryotes and eukaryotes. PPCDC (EC 4.1.1.36), one of the last enzymes in this pathway, converts phosphopantothenoylcysteine to 4-prime-phosphopantetheine (Daugherty et al., 2002 [PubMed 11923312]).
Applications:	ELISA, IHC
Name of antibody:	PPCDC
Immunogen:	Full length fusion protein
Full name:	phosphopantothenoylcysteine decarboxylase
Synonyms:	coaC; MDS018; PPC-DC

SwissProt:	Q96CD2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

