

兔抗 SLC40A1 多克隆抗体

中文名称：兔抗 SLC40A1 多克隆抗体

英文名称：Anti-SLC40A1 rabbit polyclonal antibody

别名：solute carrier family 40 member 1; FPN1; HFE4; MTP1; IREG1; MST079; MSTP079; SLC11A3

抗原：SLC40A1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a cell membrane protein that may be involved in iron export from duodenal epithelial cells. Defects in this gene are a cause of hemochromatosis type 4 (HFE4).
Applications:	ELISA, IHC
Name of antibody:	SLC40A1
Immunogen:	Synthetic peptide of human SLC40A1
Full name:	solute carrier family 40 member 1
Synonyms:	FPN1; HFE4; MTP1; IREG1; MST079; MSTP079; SLC11A3

SwissProt:	Q9NP59
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-300

