

兔抗 SLC41A3 多克隆抗体

中文名称：兔抗 SLC41A3 多克隆抗体

英文名称：Anti-SLC41A3 rabbit polyclonal antibody

别名：solute carrier family 41 member 3; SLC41A1-L2

抗原：SLC41A3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	SLC41A3 (solute carrier family 41 member 3) is a 507 amino acid multi-pass membrane protein that belongs to the SLC41A transporter family that includes SLC41A 1 and SLC41A2. Localizing to the cell membrane, SLC4 1A3 contains ten transmembrane domains and is belie ved to function as a plasma-membrane cation transpor ter. Its other two family members, SLC41A1 and SLC41 A2, have been characterized to function as magnesium transporters and are essential for maintaining magnesi um homeostasis within the cell.
Applications:	ELISA, IHC

Name of antibody:	SLC41A3
Immunogen:	Synthetic peptide of human SLC41A3
Full name:	solute carrier family 41 member 3
Synonyms:	SLC41A1-L2
SwissProt:	Q96GZ6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150

