

兔抗 PPHLN1 多克隆抗体

中文名称：兔抗 PPHLN1 多克隆抗体

英文名称： Anti-PPHLN1 rabbit polyclonal antibody

别名： HSPC206; HSPC232

抗原： PPHLN1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is one of the several proteins that become sequentially incorporated into the cornified cell envelope during the terminal differentiation of keratinocyte at the outer layers of epidermis. This protein interacts with periplakin, which is known as a precursor of the cornified cell envelope. The cellular localization pattern and insolubility of this protein suggest that it may play a role in epithelial differentiation and contribute to epidermal integrity and barrier formation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed.

Applications:	WB, IHC
Name of antibody:	PPHLN1
Immunogen:	Synthesized peptide derived from internal of human PPHLN .
Full name:	periphilin 1
Synonyms :	HSPC206; HSPC232
SwissProt:	Q8NEY8
IHC positive control:	Human liver carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	53 kDa
WB Positive control:	COLO cells and 293 cells lysates
WB Recommended dilution:	500-3000



