

## 兔抗 SLC44A2 多克隆抗体

中文名称：兔抗 SLC44A2 多克隆抗体

英文名称： Anti-SLC44A2 rabbit polyclonal antibody

别名： CTL2; PP1292

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC44A2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	CTL2, also called SLC44A2 (solute carrier family 44 member 2), is a multi-pass membrane protein expressed in cells of the inner ear. CTL2 is a possible candidate for autoimmune hearing loss in humans. Present in supporting cells of the inner ear (at protein level). Only isoform 3 is expressed in inner ear vestibular tissue.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC44A2
<b>Immunogen:</b>	Synthetic peptide of human SLC44A2
<b>Full name:</b>	solute carrier family 44 (choline transporter), member 2

<b>Synonyms :</b>	CTL2; PP1292
<b>SwissProt:</b>	Q8IWA5
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human prostate cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

