

兔抗 KRT4 多克隆抗体

中文名称: 兔抗 KRT4 多克隆抗体

英文名称: Anti-KRT4 rabbit polyclonal antibody

别名: K4, CK4, CK-4, CYK4

抗原: KRT4

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in differentiated layers of the mucosal and esophageal epithelia with family member KRT13.
Applications:	ELISA, WB
Name of antibody:	KRT4

Immunogen:	Fusion protein of human KRT4
Full name:	keratin 4, type II
Synonyms :	K4; CK4; CK-4; CYK4; WSN1
SwissProt:	P19013
ELISA Recommended dilution:	2000-10000
WB Predicted band size:	57 kDa
WB Positive control:	Hela cell
WB Recommended dilution:	1000-5000

