

兔抗 KRT6A 多克隆抗体

- 中文名称：兔抗 KRT6A 多克隆抗体
- 英文名称：Anti-KRT6A rabbit polyclonal antibody
- 别名：K6A, K6C, K6D, CK6A, CK6C, CK6D, KRT6C, KRT6D
- 抗原：KRT6A
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. As many as six of this type II cytokeratin (KRT6) have been identified; the multiplicity of the genes is attributed to successive gene duplication events. The genes are expressed with family members KRT16 and/or KRT17 in the filiform papillae of the tongue, the stratified epithelial lining of oral mucosa and esophagus, the outer root sheath of hair follicles.

	icles, and the glandular epithelia. This KRT6 gene in particular encodes the most abundant isoform. Mutations in these genes have been associated with pachyonychia congenita. The type I I cytokeratins are clustered in a region of chromosome 12q12-q13.
Applications:	ELISA, IHC
Name of antibody:	KRT6A
Immunogen:	Fusion protein of human KRT6A
Full name:	Keratin, type II cytoskeletal 6A
Synonyms :	K6A, K6C, K6D, CK6A, CK6C, CK6D, KRT6C, KRT6D
SwissProt:	P02538
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

