

兔抗 PPIG 多克隆抗体

中文名称: 兔抗 PPIG 多克隆抗体

英文名称: Anti-PPIG rabbit polyclonal antibody

别名: CYP; SRCyp; SCAF10; CARS-Cyp

抗原: PPIG

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. May be implicated in the folding, transport, and assembly of proteins. May play an important role in the regulation of pre-mRNA splicing.
Applications:	ELISA, IHC
Name of antibody:	PPIG
Immunogen:	Fusion protein of human PPIG
Full name:	peptidylprolyl isomerase G
Synonyms:	CYP; SRCyp; SCAF10; CARS-Cyp

SwissProt:	Q13427
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	35-200

