

## 兔抗 KRT8 多克隆抗体

- 中文名称：兔抗 KRT8 多克隆抗体
- 英文名称：Anti-KRT8 rabbit polyclonal antibody
- 别名：K8; KO; CK8; CK-8; CYK8; K2C8; CARD2
- 储存：冷冻（-20℃）
- 抗原：KRT8
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene is a member of the type II keratin family clustered on the long arm of chromosome 12. Type I and type II keratins heteropolymerize to form intermediate-sized filaments in the cytoplasm of epithelial cells. The product of this gene typically dimerizes with keratin 18 to form an intermediate filament in simple single-layered epithelial cells. This protein plays a role in maintaining cellular structural integrity and also functions in signal transduction and cellular differentiation. Mutations in this gene cause cryptogenic cirrhosis. Alternatively spliced transcript variants have been found for

	this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KRT8
<b>Immunogen:</b>	Fusion protein of human KRT8
<b>Full name:</b>	keratin 8
<b>Synonyms:</b>	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2
<b>SwissProt:</b>	P05787
<b>IHC positive control:</b>	Human gastric cancer and human breast cancer
<b>IHC Recommend dilution:</b>	25-100

