

## 兔抗 ZNF34 多克隆抗体

中文名称：兔抗 ZNF34 多克隆抗体

英文名称：Anti-ZNF34 rabbit polyclonal antibody

别名：KOX32

抗原：ZNF34

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ZNF34
<b>Immunogen:</b>	Fusion protein of human ZNF34
<b>Full name:</b>	zinc finger protein 34
<b>Synonyms :</b>	KOX32
<b>SwissProt:</b>	Q8IZ26
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer and human breast cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	64 kDa

<b>WB Positive control:</b>	Human liver cancer tissue lysate
<b>WB Recommended dilution:</b>	200-1000



