

## 兔抗 KRT9 多克隆抗体

中文名称：兔抗 KRT9 多克隆抗体

英文名称：Anti-KRT9 rabbit polyclonal antibody

别名：keratin 9 K9; CK-9; EPPK

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KRT9

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the type I keratin 9, an intermediate filament chain expressed only in the terminally differentiated epidermis of palms and soles. Mutations in this gene cause epidermolytic palmoplantar keratoderma.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KRT9
<b>Immunogen:</b>	Synthetic peptide of human KRT9
<b>Full name:</b>	keratin 9

<b>Synonyms:</b>	K9; CK-9; EPPK
<b>SwissProt:</b>	P35527
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	30-150



