

兔抗 KRT12 多克隆抗体

中文名称: 兔抗 KRT12 多克隆抗体

英文名称: Anti-KRT12 rabbit polyclonal antibody

别名: keratin 12; K12

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KRT12

反应种属: Human, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	KRT12 encodes the type I intermediate filament chain keratin 12, expressed in corneal epithelia. Mutations in this gene lead to Meesmann corneal dystrophy.
Applications:	ELISA, WB
Name of antibody:	KRT12
Immunogen:	Synthetic peptide of human KRT12
Full name:	keratin 12
Synonyms:	K12
SwissProt:	Q99456

ELISA Recommended dilution:	5000-10000
WB Predicted band size:	54 kDa
WB Positive control:	Rat eye tissue and Human skin tissue lysates
WB Recommended dilution:	200-1000

