

兔抗 SLITRK4 多克隆抗体

中文名称：兔抗 SLITRK4 多克隆抗体

英文名称： Anti-SLITRK4 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SLITRK4

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane protein belonging to the the SLITRK family. These family members include two N-terminal leucine-rich repeat domains similar to those found in the axonal growth-controlling protein SLIT, as well as C-terminal regions similar to neurotrophin receptors. Studies of an homologous protein in mouse suggest that this family member functions to suppress neurite outgrowth. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	SLITRK4
Immunogen:	Fusion protein of human SLITRK4
Full name:	SLIT and NTRK-like family, member 4

SwissProt:	Q8IW52
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	10-50
WB Predicted band size:	94 kDa
WB Positive control:	PC3 cells
WB Recommended dilution:	200-1000



