

## 兔抗 KRT13 多克隆抗体

中文名称：兔抗 KRT13 多克隆抗体

英文名称：Anti-KRT13 rabbit polyclonal antibody

别名：K13; CK13

相关类别：一抗

储存：冷冻（-20℃）

抗原：KRT13

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the suprabasal layers of non-cornified stratified epithelia. Mutations in this gene and keratin 4 have been associated with the autosomal dominant disorder White Sponge Nevus. The typ

	<p>e I cytokeratins are clustered in a region of chromosome 17 q21.2. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been described.</p>
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KRT13
<b>Immunogen:</b>	Fusion protein of human KRT13
<b>Full name:</b>	keratin 13
<b>Synonyms :</b>	K13; CK13
<b>SwissProt:</b>	P13646
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human lung cancer
<b>IHC Recommend dilution:</b>	50-200

