

兔抗 SLITRK5 多克隆抗体

中文名称：兔抗 SLITRK5 多克隆抗体

英文名称：Anti-SLITRK5 rabbit polyclonal antibody

别名：LRRC11; bA364G4.2

抗原：SLITRK5

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Members of the SLITRK family, such as SLITRK5, are integral membrane proteins with 2 N-terminal leucine-rich repeat (LRR) domains similar to those of SLIT proteins (see SLIT1; MIM 603742). Most SLITRKs, including SLITRK5, also have C-terminal regions that share homology with neurotrophin receptors (see NTRK1; MIM 191315). SLITRKs are expressed predominantly in neural tissues and have neurite-modulating activity (Aruga et al., 2003 [PubMed 14557068]).
Applications:	ELISA, IHC

Name of antibody:	SLITRK5
Immunogen:	Fusion protein of human SLITRK5
Full name:	SLIT and NTRK like family member 5
Synonyms:	LRRC11; bA364G4.2
SwissProt:	O94991
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human prostate cancer
IHC Recommend dilution:	25-100



