

兔抗 KRT15 多克隆抗体

中文名称: 兔抗 KRT15 多克隆抗体

英文名称: Anti-KRT15 rabbit polyclonal antibody

别名: K15, CK15, K1CO

储存: 冷冻 (-20℃)

抗原: KRT15

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains and are clustered in a region on chromosome 17q21.2.
Applications:	ELISA, WB
Name of antibody:	KRT15
Immunogen:	Fusion protein of human KRT15

Full name:	Keratin, type I cytoskeletal 15
Synonyms :	K15, CK15, K1CO
SwissProt:	P19012
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	49 kDa
WB Positive control:	Human ovarian cancer tissue
WB Recommended dilution:	500-2000

