

## 兔抗 SMAD1 (Phospho-Ser465)多克隆抗体

中文名称：兔抗 SMAD1 (Phospho-Ser465)多克隆抗体

英文名称：Anti-SMAD1 (Phospho-Ser465) rabbit polyclonal antibody

别名：BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SMAD1 (Phospho-Ser465)

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcriptional modulator activated by BMP (bone morphogenetic proteins) type 1 receptor kinase. SMAD1 is a receptor-regulated SMAD (R-SMAD).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SMAD1 (Phospho-Ser465)
<b>Immunogen:</b>	Synthetic peptide of human SMAD1 (Phospho-Ser465)
<b>Full name:</b>	SMAD family member 1 (Phospho-Ser465)
<b>Synonyms :</b>	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1
<b>SwissProt:</b>	Q15797
<b>WB Predicted band size:</b>	52 kDa

<b>WB Positive control:</b>	Hela cells treated with PMA
<b>WB Recommended dilution:</b>	500-1000

