

## 兔抗 KRT17 多克隆抗体

中文名称：兔抗 KRT17 多克隆抗体

英文名称：Anti-KRT17 rabbit polyclonal antibody

别名：PC, K17, PC2, PCHC1

储存：冷冻（-20℃）

抗原：KRT17

宿主：Rabbit

反应种属 Human, Mouse, Rat

相关类别 一抗

标记物 Unconjugate

克隆类型 rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	KRT17
<b>Immunogen:</b>	Fusion protein of human KRT17
<b>Full name:</b>	keratin 17
<b>Synonyms:</b>	PC; K17; PC2; 39.1; CK-17; PCHC1

<b>SwissProt:</b>	Q04695
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	48 kDa
<b>WB Positive control:</b>	SKOV3 and HEPG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000





