

## 兔抗 SMAD1(Ab-187) 多克隆抗体

中文名称: 兔抗 SMAD1(Ab-187) 多克隆抗体

英文名称: Anti-SMAD1(Ab-187) rabbit polyclonal antibody

别 名: SMAD family member 1

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: SMAD1(Ab-187)

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

**Background:** 

The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein mediates the signals of the bone morphogenetic proteins (BMPs), which are involved in a range of biological activities including cell growth, apoptosis, morphogenesis, development and immune respon



	ses.
Applications:	WB, IF
Name of antibody:	SMAD1(Ab-187)
Immunogen:	Synthesized non-phosphopeptide derived from h uman Smad1 around the phosphorylation site of serine 187 (P-H-S(p)-P-N).
Full name:	SMAD family member 1
Synonyms:	SMAD family member 1
SwissProt:	Q15797
WB Predicted band size:	52 kDa
WB Positive control:	Jurkatlysates K562 cell
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500





