

兔抗 SMAD1(Ab-465) 多克隆抗体

中文名称: 兔抗 SMAD1(Ab-465) 多克隆抗体

英文名称: Anti-SMAD1(Ab-465) rabbit polyclonal antibody

别名: BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: SMAD1(Ab-465)

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|----------------------|--|
| Background: | The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein mediates the signals of the bone morphogenetic proteins (BMPs), which are involved in a range of biological activities including cell growth, apoptosis, morphogenesis, development and immune responses. |
| Applications: | WB, IHC |

| | |
|---------------------------------|---|
| Name of antibody: | SMAD1(Ab-465) |
| Immunogen: | Synthesized non-phosphopeptide derived from human Smad1 around the phosphorylation site of serine 465 (I-S-S-V-S(p)). |
| Full name: | SMAD family member 1 |
| Synonyms : | BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1 |
| SwissProt: | Q15797 |
| IHC positive control: | Human breast carcinoma tissue |
| IHC Recommend dilution: | 50-100 |
| WB Predicted band size: | 52 kDa |
| WB Positive control: | COS7 cells treated with EGF (200ng/ml, 15mins) lysate |
| WB Recommended dilution: | 500-3000 |



