

## 兔抗 PPP1R14C 多克隆抗体

中文名称：兔抗 PPP1R14C 多克隆抗体

英文名称： Anti-PPP1R14C rabbit polyclonal antibody

别名： KEPI; NY-BR-81; CPI17-like

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PPP1R14C

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The degree of protein phosphorylation is regulated by a balance of protein kinase and phosphatase activities. Protein phosphatase-1 (PP1; see MIM 176875) is a signal-transducing phosphatase that influences neuronal activity, protein synthesis, metabolism, muscle contraction, and cell division. PPP1R14C is an inhibitor of PP1 (Liu et al., 2002 [PubMed 11812771]).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PPP1R14C

<b>Immunogen:</b>	Synthesized peptide derived from internal of human PP14 C.
<b>Full name:</b>	protein phosphatase 1, regulatory (inhibitor) subunit 14C
<b>Synonyms :</b>	KEPI; NY-BR-81; CPI17-like
<b>SwissProt:</b>	Q8TAE6
<b>WB Predicted band size:</b>	18 kDa
<b>WB Positive control:</b>	Jurkat cells and COLO205 lysates
<b>WB Recommended dilution:</b>	500-3000

