

兔抗 PPP1R15A 多克隆抗体

中文名称：兔抗 PPP1R15A 多克隆抗体

英文名称： Anti-PPP1R15A rabbit polyclonal antibody

别名： GADD34

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PPP1R15A

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The induction of this gene by ionizing radiation occurs in certain cell lines regardless of p53 status, and its protein response is correlated with apoptosis following ionizing radiation.
Applications:	ELISA, IHC
Name of antibody:	PPP1R15A
Immunogen:	Synthetic peptide of human PPP1R15A
Full name:	protein phosphatase 1, regulatory subunit 15A

Synonyms :	GADD34
SwissProt:	O75807
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	30-150

