

兔抗 PPP1R16A 多克隆抗体

中文名称: 兔抗 PPP1R16A 多克隆抗体

英文名称: Anti-PPP1R16A rabbit polyclonal antibody

别名: MYPT3

抗原: PPP1R16A

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Inhibits protein phosphatase 1 activity toward phosphorylase , myosin light chain and myosin substrates By similarity.
Applications:	WB
Name of antibody:	PPP1R16A
Immunogen:	Synthesized peptide derived from internal of human PPP1R16A.
Full name:	protein phosphatase 1, regulatory subunit 16A
Synonyms :	MYPT3
SwissProt:	Q96I34

WB Predicted band size:	58 kDa
WB Positive control:	Rat heart tissue lysate
WB Recommended dilution:	500-3000

