

## 兔抗 PPP1R1A 多克隆抗体

中文名称：兔抗 PPP1R1A 多克隆抗体

英文名称： Anti-PPP1R1A rabbit polyclonal antibody

别名： protein phosphatase 1 regulatory inhibitor subunit 1A; I1; IPP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PPP1R1A

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Inhibitor of protein-phosphatase 1. This protein may be important in hormonal control of glycogen metabolism. Hormones that elevate intracellular cAMP increase I-1 activity in many tissues. I-1 activation may impose cAMP control over proteins that are not directly phosphorylated by PKA. Following a rise in intracellular calcium, I-1 is inactivated by calcineurin (or PP2B). Does not inhibit type-2 phosphatases.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PPP1R1A

<b>Immunogen:</b>	Fusion protein of human PPP1R1A
<b>Full name:</b>	protein phosphatase 1 regulatory inhibitor subunit 1A
<b>Synonyms:</b>	I1; IPP1
<b>SwissProt:</b>	Q13522
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human brain
<b>IHC Recommend dilution:</b>	50-300



