

## 兔抗 PPP1R1B 多克隆抗体

中文名称：兔抗 PPP1R1B 多克隆抗体

英文名称： Anti-PPP1R1B rabbit polyclonal antibody

别名： DARPP32; DARPP-32

相关类别： 一抗

抗原： PPP1R1B

储存： 冷冻 (-20℃)

宿主： Rabbit

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a bifunctional signal transduction molecule. Dopaminergic and glutamatergic receptor stimulation regulates its phosphorylation and function as a kinase or phosphatase inhibitor. As a target for dopamine, this gene may serve as a therapeutic target for neurologic and psychiatric disorders. Multiple transcript variants encoding different isoforms have been found for this gene
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PPP1R1B

<b>Immunogen:</b>	Fusion protein of human PPP1R1B
<b>Full name:</b>	protein phosphatase 1, regulatory (inhibitor) subunit 1B
<b>Synonyms :</b>	DARPP32; DARPP-32
<b>SwissProt:</b>	Q9UD71
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human colon cancer
<b>IHC Recommend dilution:</b>	50-200

