

兔抗 PPP1R42 多克隆抗体

中文名称: 兔抗 PPP1R42 多克隆抗体

英文名称: Anti-PPP1R42 rabbit polyclonal antibody

别名: dtr; TLLR; LRRC67

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PPP1R42

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Regulates phosphatase activity of protein phosphatase 1 (PP1) complexes in the testis.
Applications:	ELISA, IHC
Name of antibody:	PPP1R42
Immunogen:	Synthetic peptide of human PPP1R42
Full name:	protein phosphatase 1 regulatory subunit 42
Synonyms:	dtr; TLLR; LRRC67
SwissProt:	Q7Z4L9
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100



