

兔抗 PPP1R7 多克隆抗体

中文名称: 兔抗 PPP1R7 多克隆抗体

英文名称: Anti-PPP1R7 rabbit polyclonal antibody

别名: SDS22

抗原: PPP1R7

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a protein subunit that regulates the activity of the serine/threonine phosphatase, protein phosphatase-1. The encoded protein is required for completion of the mitotic cycle and for targeting protein phosphatase-1 to mitotic kinetochores. Alternate splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	PPP1R7
Immunogen:	Synthesized peptide derived from internal of human PPP1R7.

Full name:	protein phosphatase 1, regulatory subunit 7
Synonyms :	SDS22
SwissProt:	Q15435
WB Predicted band size:	42 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000

