

兔抗 PPP2CB 多克隆抗体

中文名称: 兔抗 PPP2CB 多克隆抗体

英文名称: Anti-PPP2CB rabbit polyclonal antibody

别名: PP2CB, PP2Abeta

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PPP2CB

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes a beta isoform of the catalytic subunit.
Applications:	ELISA, IHC
Name of antibody:	PPP2CB

Immunogen:	Synthetic peptide of human PPP2CB
Full name:	Protein phosphatase 2, catalytic subunit, beta isozyme
Synonyms :	PP2CB, PP2Abeta
SwissProt:	P62714
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and Human tonsil
IHC Recommend dilution:	100-300

