

兔抗 PPP2R3B 多克隆抗体

中文名称：兔抗 PPP2R3B 多克隆抗体

英文名称： Anti-PPP2R3B rabbit polyclonal antibody

别名： PR48; NYREN8; PPP2R3L; PPP2R3LY

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PPP2R3B

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

Protein phosphatase 2 (formerly named type 2A) is one of the four major Ser/Thr phosphatases and is implicated in the negative control of cell growth and division. Protein phosphatase 2 holoenzymes are heterotrimeric proteins composed of a structural subunit A, a catalytic subunit C, and a regulatory subunit B. The regulatory subunit is encoded by a diverse set of genes that have been grouped into the B/PR55, B'/PR61, and B''/PR72 families. These different regulatory subunits confer distinct enzymatic specificities and intracellular localizations to the holoenzyme. The product of

	this gene belongs to the B'' family. The B'' family has been further divided into subfamilies. The product of this gene belongs to the beta subfamily of regulatory subunit B''.
Applications:	WB
Name of antibody:	PPP2R3B
Immunogen:	Synthesized peptide derived from internal of human P2R3B .
Full name:	protein phosphatase 2, regulatory subunit B'', beta
Synonyms :	PR48; NYREN8; PPP2R3L; PPP2R3LY
SwissProt:	Q9Y5P8
WB Predicted band size:	65 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

