

## 兔抗 PPRC1 多克隆抗体

中文名称: 兔抗 PPRC1 多克隆抗体

英文名称: Anti-PPRC1 rabbit polyclonal antibody

别名: PRC

抗原: PPRC1

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is similar to PPAR-gamma coactivator 1 (PPARGC1/PGC-1), a protein that can activate mitochondrial biogenesis in part through a direct interaction with nuclear respiratory factor 1 (NRF1). This protein has been shown to interact with NRF1. It is thought to be a functional relative of PPAR-gamma coactivator 1 that activates mitochondrial biogenesis through NRF1 in response to proliferative signals. Alternative splicing results in multiple transcript variants
<b>Applications:</b>	IHC

<b>Name of antibody:</b>	PPRC1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human PPRC1.
<b>Full name:</b>	peroxisome proliferator-activated receptor gamma, coactivator-related 1
<b>Synonyms :</b>	PRC
<b>SwissProt:</b>	Q5VV67
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100

