

## 兔抗 PRELID1 多克隆抗体

中文名称：兔抗 PRELID1 多克隆抗体

英文名称：Anti-PRELID1 rabbit polyclonal antibody

别名：PX19; PRELI; SBBI12; CGI-106

相关类别：一抗

抗原：PRELID1

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the late embryogenesis abundant motif-containing protein family. The encoded protein is localized to mitochondria and may function as a cytoprotectant by regulating cell death and differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. Several related pseudogenes have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PRELID1
<b>Immunogen:</b>	Fusion protein of human PRELID1
<b>Full name:</b>	PRELI domain containing 1
<b>Synonyms :</b>	PX19; PRELI; SBBI12; CGI-106

<b>SwissProt:</b>	Q9Y255
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	30-150

