

## 兔抗 PREPL 多克隆抗体

中文名称：兔抗 PREPL 多克隆抗体

英文名称： Anti-PREPL rabbit polyclonal antibody

别名： prolyl endopeptidase like; CMS22

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PREPL

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PREPL
<b>Immunogen:</b>	Fusion protein of human PREPL

<b>Full name:</b>	prolyl endopeptidase like
<b>Synonyms:</b>	CMS22
<b>SwissProt:</b>	Q4J6C6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	84 kDa
<b>WB Positive control:</b>	RAW264.7 cell lysate
<b>WB Recommended dilution:</b>	500-2000



