

## 兔抗 PRKCA(Phospho-Tyr657) 多克隆抗体

中文名称: 兔抗 PRKCA(Phospho-Tyr657) 多克隆抗体

英文名称: Anti-PRKCA(Phospho-Tyr657) rabbit polyclonal antibody

别 名: AAG6; PKCA; PRKACA; PKC-alpha

抗原: PRKCA(Phospho-Tyr657)

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

**Background:** 

Protein kinase C (PKC) is a family of serine- and threonine -specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family mem bers phosphorylate a wide variety of protein targets and a re known to be involved in diverse cellular signaling path ways. PKC family members also serve as major receptors f or phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role in cells. The protein encoded by this gene is one of the PKC family members.



Applications:	WB
Name of antibody:	PRKCA(Phospho-Tyr657)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 657(F-S-Y(p)-V-N) derived from Human PKC $\alpha$ .
Full name:	protein kinase C, alpha
Synonyms:	AAG6; PKCA; PRKACA; PKC-alpha
SwissProt:	P17252
WB Predicted band size:	77 kDa
WB Positive control:	COLO205 lysates
WB Recommended dilution:	500-1000

